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Table 2. HiTrap Desalting characteristics

Matrix	Sephadex G-25 Superfine, cross-linked dextran
Separation mechanism	According to size
Void volume	1.5 ml
Recommended sample volume	0.1 to 1.5 ml
Sample dilution, syringe operation	1.3 to 4.0 × applied volume
Exclusion limit	M _r 5 000
Bead size	15 to 70 µm
Maximum flow rate ¹	15 ml/min
Recommended flow rate ¹	1 to 10 ml/min
Back pressure at 10 ml/min ¹	0.25 bar
Chemical stability	All commonly used buffers
pH stability, short and long term ²	2 to 13
Storage	20% ethanol at 4°C to 30°C

¹ room temperature, aqueous buffers

² *short term* refers to the pH interval for regeneration.

long term refers to the pH interval where the medium is stable over a long period of time without adverse effects on its subsequent chromatographic performance.

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Table 3. Recommended sample and elution volumes using a syringe. Examples of typical yields and remaining salt in the desalted sample.

Sample load (ml)	Elute and Add buffer (ml)	Collect (ml)	Remaining Yield (%)	Salt (%)	Dilution factor
0.25	1.25	1.0	> 95	0.0	4.0
0.50	1.0	1.5	> 95	< 0.1	3.0
1.00	0.5	2.0	> 95	< 0.2	2.0
1.50	0.0	2.0	> 95	< 0.2	1.3

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Column: HiTrap Desalting
Sample: 2 mg/ml bovine serum albumin in 50 mM sodium phosphate buffer, 0.5 M sodium chloride, pH 7.0
Sample volume: 0.3 ml, 1.3 ml, 2.2 ml
Eluent: 50 mM sodium phosphate buffer, 0.15 M sodium chloride, pH 7.0
Flow rate: 5 ml/min
Detection: UV (280 nm, 5 mm cell) and conductivity
Sample Injection: Syringe


Fig 2. The effect of different sample volumes on the HiTrap Desalting column using a syringe for sample injection

Column:	HiTrap Desalting
Sample:	2 mg/ml bovine serum albumin in 50 mM sodium phosphate buffer, 0.5 M sodium chloride, pH 7.0
Sample volume:	0.5 ml, 1.1 ml, 2.1 ml
Eluent:	50 mM sodium phosphate buffer, 0.15 M sodium chloride, pH 7.0
Flow rate:	5 ml/min
Detection:	UV (280 nm, 5 mm cell) and conductivity
Sample injection:	Tubing loops


Fig 3. The effect of different sample volumes on the HiTrap Desalting column using tubing loops for sample injection.

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Accessories	Quantity	Code No.
1/16" male/luer female <i>(For connection of syringe to top of HiTrap column)</i>	2	18-1112-51
Tubing connector flangeless/M6 female <i>(For connection of tubing to bottom of HiTrap column)</i>	2	18-1003-68
Tubing connector flangeless/M6 male <i>(For connection of tubing to top of HiTrap column)</i>	2	18-1017-98
Union 1/16" female/M6 male <i>(For connection to original FPLC System through bottom of HiTrap column)</i>	6	18-1112-57
Union M6 female /1/16" male <i>(For connection to original FPLC System through top of HiTrap column)</i>	5	18-3858-01
Union luerlock female/M6 female	2	18-1027-12
HiTrap/HiPrep, 1/16" male connector for ÄKTA design	8	28-4010-81
Stop plug female, 1/16" <i>(For sealing bottom of HiTrap column)</i>	5	11-0004-64
Fingertight stop plug, 1/16"	5	11-0003-55